

ZA070A Zero Air

発電機

品番: UK: 30-1610 (230v) | EU: 30-1620 (230v) | US: 30-1630 (110v)



説明

ZA070Aジェネレーターは、GCアプリケーションで検出器ガスのゼロエアを生成し、機器に入る前に圧縮空気源からすべての炭化水素を除去して、検出器からのきれいな結果を保証します。




Generator Details

- ガスの種類: 空気
- 統合型コンプレッサ: 番号

主な機能

- 0.1 ppm未満の炭化水素含有量で、ゼログレードエアの一定流量を生成します
- 必要に応じて24時間年中無休で実施
- 1つの発電機で複数の機器に供給できる
- 可動部品が少ないため、メンテナンスがほとんど不要で、発電機の寿命が長くなります
- 最小限の設定が必要
- 12か月の包括的なオンサイト保証



技術仕様

ZA070A Zero Air 発電機	
ガスの種類	空気
最大ガス流量:	7L/min
炭化水素の含有量:	0.1ppm
最高純度:	0.1ppm THC
起動時間:	40 Minutes
電圧:	230 ± 10% VAC
周波数:	50 HZ
電流:	4.9 Amps
動作時の最高温度:	30°C / 86°F
分子:	0.01µm
大きさ (高さx幅x奥行き) ミリメートル	880 x 430 x 410 mm
大きさ (高さx幅x奥行き) インチ	34.6 x 16.9 x 16.1 inches
装置の重量	64kg / 140.8lbs

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001

