

Precision Zero Air 3.5L


ガスジェネレータ

品番: 60-3500 (230V) | 60-3501 (110V)



説明

Precision Zero Air発生装置は、標準検出限界と微量検出限界レベルの両方で、GCのフレイムサポートガスとして使用される炭化水素フリーの清潔で乾燥したエアを供給するように特別に設計されています




Generator Details

- ガスの種類: 空気
- 統合型コンプレッサ: 番号

主な機能

- 他のPrecisionユニットとシームレスに統合します
- 炭化水素含有量(メタンとして) <math>< 0.05\text{ppm}</math>市場をリードする純度
- 最小限の寿命維持要件、高価な触媒チャンバーの交換なし
- システムに汚染物質が入るリスクを回避する(空のシリンダーを交換する場合)
- 分析の途中でガスが不足するリスクがない
- Precision Airに対応
- コンプレッサーの状態を示す照明機能(黄色のウォームアップ、緑の準備完了)
- 12か月の包括的なオンサイト保証



技術仕様

	Precision Zero Air 3.5L - ガスジェネレータ
ガスの種類	空気
最大ガス流量:	3.5L/min
最大出口圧力:	80psi/5.5bar
炭化水素の含有量:	<0.05ppm
最大出口圧力:	80psi/5.5bar
ガス出口の設定:	1 x 1/4" BSPP
最低入口圧力:	90psi/6.2bar
最高入口圧力:	145psi/10bar
Min Air Inlet Flow:	3.5 L/min
最小入口空気品質:	ISO 8753-1:2010 Class 1.4.1 or Precision Air Compressor
起動時間:	60 Minutes
電圧:	110 / 230 ± 10% V
周波数:	50 / 60 Hz
電流:	1 - 2 A
熱出力:	785
動作時の最高温度:	35°C / 95°F
認証:	CE, FCC, CSA
大きさ (高さx幅x奥行き) ミリメートル	156 x 380 x 540 mm
大きさ (高さx幅x奥行き) インチ	6.1 x 15 x 21.3 inches
装置の重量	16kg / 35.2lbs

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001

