

Precision Air Compressor


精密レンジ用エアコンプレッサー

品番: 65-0555 (230V) | 65-1555 (110V)



説明

Precision Air Compressorは、Precision NitrogenおよびPrecision Zero Airガスジェネレーターに専用の完全に統合された圧縮空気源を提供するように特別に設計および設計されています。発電機と同じコンパクトで積み重ね可能なフォームファクター内に収容されたPrecision Air Compressorモジュールは、通常はベースユニットとして任意のPrecisionスタックに追加できます。利用可能です。独自の二重コンプレッサー取り付けシステムのおかげで、Precision Air Compressorは低ノイズと低振動のメリットも享受し、通常の実験室での操作において比較的静かになります。




Generator Details

- ガスの種類: 空気
- 統合型コンプレッサ: はい

主な機能

- 精密窒素およびゼロエアジェネレーター用の圧縮空気供給
- 1つのユニットで、Precision NitrogenとZero Airの両方の発電機をスタックで供給できます
- 異音や振動を減衰させる二重コンプレッサー取り付けシステム
- サービス&修理ステータスインジケータ
- 12か月の包括的なオンサイト保証



技術仕様

	Precision Air Compressor
ガスの種類	空気
最大ガス流量:	VariousL/min
最大出口圧力:	120psi/8.3bar
最大出口圧力:	120psi/8.3bar
ガス出口の設定:	1 x 1/4" BSPP
消費電力:	510 / 1150 W
電圧:	110 / 230 V
周波数:	50 / 60 Hz
電流:	10 A (110V) / 2.45 A (230V)
熱出力:	1740 - 3922 (BTU)
動作時の最高温度:	35°C / 95°F
認証:	CSA, FCC, CE
大きさ (高さx幅x奥行き) ミリメートル	406 x 380 x 540 mm
大きさ (高さx幅x奥行き) インチ	16 x 15 x 21.3 inches
装置の重量	42kg / 92.4lbs

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001

