



## 技術仕様

	Precision Zero Air 30 02発電機
ガスの種類	空気
最大ガス流量:	30L/min
最大出口圧力:	100psi/6.9bar
炭化水素の含有量:	<0.05ppm
最大出口圧力:	100psi/6.9bar
ガス出口の設定:	1 x 1/4" BSPP
最低入口圧力:	110psi/7.6bar
最高入口圧力:	145psi/10bar
Min Air Inlet Flow:	30 L/min
最小入口空気品質:	ISO 8573-1:2010 Class 1.4.1 or Precision Air Compressor
起動時間:	60 mins
消費電力:	1210/1541 W
電圧:	110 / 230 V
周波数:	50 / 60 Hz
電流:	6.7 - 11 A
熱出力:	5255
動作時の最高温度:	35°C / 95°F
認証:	CE, FCC, CSA
大きさ (高さx幅x奥行き) ミリメートル	406 x 380 x 540 mm
大きさ (高さx幅x奥行き) インチ	16 x 15 x 21.3 inches
装置の重量	41kg / 90.2lbs

# [PEAK Protected]<sup>TM</sup>

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: [www.peakscientific.com/ordering/](http://www.peakscientific.com/ordering/)
- For service plans visit: [www.peakscientific.com/service/service-plans/](http://www.peakscientific.com/service/service-plans/)

**Peak Scientific UK**  
Tel: +44 (0)141 812 8100  
Fax: +44 (0)141 812 8200

**Peak Scientific North America**  
Tel: +1-800-767-6532  
Fax: +1-978-608-9503

**Peak Scientific China**  
Tel: +86 21 5079 1190  
Fax: +86 21 5079 1191  
For a full list of our worldwide office locations, please visit:  
Web: [www.peakscientific.com](http://www.peakscientific.com)  
Email: [discover@peakscientific.com](mailto:discover@peakscientific.com)

Peak Scientific's Quality Management System conforms to: ISO:9001

