

TOC 1500

ゼロエア CO2 フリー

品番: 42-1030



説明

TOCアナライザーに信頼性の高い空気供給を提供するように設計されたTOC1500HPは、ラボでクリーンで乾燥したCO2およびCH4を生成しない堅牢なシステムです。




Generator Details

- ガスの種類: 空気
- 統合型コンプレッサ: 番号

主な機能

- クリーンで乾燥した、炭化水素と二酸化炭素を含まない空気を生成する
- TC-ICプロセスTOCアナライザーに最適
- 必要に応じて24時間年中無休で実施
- 可動部品が少ないため、メンテナンスがほとんど不要で、発電機の寿命が長くなります
- 完全に無音の動作
- 最小限の設定が必要
- Precision Air Compressorでサポート可能
- 12か月の包括的なオンサイト保証



技術仕様

	TOC 1500 Zero Air Co2 Free
ガスの種類	空気
最大ガス流量:	1500cc/min
最大出口圧力:	80psi/5.5bar
最大出口圧力:	80psi/5.5bar
最高純度:	<1 ppm CO2/ 2.0ppm THC
最低入口圧力:	90psi/6.2bar
最高入口圧力:	145psi/10bar
水質要件:	N/A
水消費量:	TBC
電圧:	100-230 ± 10% VAC
周波数:	50 / 60 Hz
電流:	8 Amps
認証:	CE
大きさ (高さx幅x奥行き) ミリメートル	660 x 250 x 160 mm
大きさ (高さx幅x奥行き) インチ	26 x 9.8 x 6.3 inches
装置の重量	12kg / 26.4lbs

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001

